





## Trust Socket

Unique transtibial diagnostic socket cone can be Direct-Fit on the patient's limb to optimize socket fit and comfort.

Low-temp, proprietary thermoplastic socket cone can be rapidly adjusted on the spot

Trust Socket has no distal adapter included so clinician can optimize alignment.

Dynamic transtibial Diagnostic socket fitting can be completed within 1-Hour!



## ...with Full Weight Bearing technology.

- Eliminates PLASTER!
- Creates a weight bearing, Direct-Fit transtibial Diagnostic socket within 45 minutes!
- Reduces the number patient visits needed to get to aligned, zero-ply Diagnostic socket fit